

## HEINE® mini-c Clip Lamp

## Compact, long-life pocket diagnostic light



Sturdy, long-life pocket diagnostic light in a compact, modern design with chrome head. Bright, concentrated light thanks to novel XHL Xenon Halogen lens-bulb.

mini-c Clip Lamp only fits mini-c handles. The mini-c handle is not compatible with mini3000 instruments. Handle available in black or blue.

- ∴ **Clip Lamp with Ear Light for otoplasty.** The Ear Light is used to ensure that hearing aids fit correctly.
- ∴ **Compact, oval, attractive light.** Lies ergonomically in the hand.
- ∴ **High-quality handle: Chrome-finish upper section/refined plastic.** Shockproof, sturdy, non-slip.
- ∴ **Attachment clip with integrated on/off switch.** Secure. Switches off automatically when replaced in the pocket. 20.000 switch cycles guaranteed.
- ∴ **Replaceable batteries.** Size AAA.
- ∴ **Registered Design.**

|   |                     |
|---|---------------------|
| <b>mini-c Clip Lamp</b>                                       | <b>2.5V</b>         |
| with batteries (in blister pack) [01]                         | <b>D-001.73.109</b> |
| with Ear Light with batteries (in blister pack) [02]          | <b>D-001.73.103</b> |
| Ear Light for mini-c Clip Lamp, 6 pieces in blister pack [03] | <b>D-000.73.105</b> |
| XHL Xenon Halogen spare bulb                                  | <b>X-001.88.108</b> |



HEINE mini3000 pocket instruments – a complete range of quality diagnostic instruments for your pocket. Available in black or blue. HEINE mini3000 Clip Lamp and mini-c Clip Lamp attachments are only available with a chrome head. Please specify the colour you prefer when ordering. If not otherwise specified, we will supply black instruments.

## HEINE mini3000® Clip Lamp

### Compact, long-life pocket diagnostic light

Sturdy, long-life pocket diagnostic lights in a compact, modern design with chrome head. Bright, concentrated light thanks to novel XHL Xenon Halogen lens-bulb.

Can only be used with the mini3000 handle system. Handle available in black or blue.



[01]

[02]

- ∴ Long-life pocket light
- ∴ Switches off automatically
- ∴ Optional rechargeable handle

- ∴ **Compact, oval, attractive light.** Lies ergonomically in the hand.
- ∴ **2-piece handle and head.** Easy to maintain, flexible. Compatible with other mini3000 instrument heads.
- ∴ **Enhanced XHL Xenon Halogen lens-bulb with 100% more light compared with conventional bulbs.** Bright, focused light.
- ∴ **High-quality handle: Chrome-finish upper section/refined plastic.** Shockproof, sturdy, non-slip.
- ∴ **Attachment clip with integrated on/off switch.** Secure. Switches off automatically when replaced in the pocket. 20.000 switch cycles guaranteed.
- ∴ **Replaceable batteries.** Size AA, or optional rechargeable battery with mini NT table charger\*.
- ∴ **Ear Light for otoplasty [02].** The Ear Light is used to ensure that hearing aids fit correctly.

| mini3000 Clip Lamp  | 2.5V         |
|---|--------------|
| complete with mini3000 battery handle, with batteries [01]      | D-001.73.131 |
| mini3000 Clip Lamp head <b>without handle</b>                   | D-001.73.130 |
| Ear Light for mini3000 Clip Lamp, 6 pieces in blister pack [02] | D-000.73.105 |
| XHL Xenon Halogen spare bulb                                    | X-001.88.107 |



\*To upgrade to a rechargeable instrument with mini NT table charger, please see page 146.

HEINE mini3000 pocket instruments – a complete range of quality diagnostic instruments for your pocket. Available in black or blue. HEINE mini3000 Clip Lamp and mini-c Clip Lamp attachments are only available with a chrome head. Please specify the colour you prefer when ordering. If not otherwise specified, we will supply black instruments.



## HEINE mini3000® Combi Lamp

### Compact, all purpose illuminator with integrated tongue-blade holder

Sturdy, long-life pocket diagnostic light in a compact, modern design. Bright, concentrated light thanks to novel XHL Xenon Halogen lens-bulb.

Can only be used with the mini3000 handle system.  
Available in black or blue.

- ∴ **Many applications.** Can be used as a general illuminator and as an illuminated tongue-blade holder for wooden blades and HEINE disposable blades.
- ∴ **2-piece handle and head.** Easy to maintain, flexible. Compatible with other mini3000 instrument heads.
- ∴ **Enhanced XHL Xenon Halogen lens-bulb with 100% more light compared with conventional bulbs.** Bright, focused light.
- ∴ **High-quality handle: Chrome-finish upper section/refined plastic.** Shockproof, sturdy, non-slip.
- ∴ **Attachment clip with integrated on/off switch.** Secure. Switches off automatically when replaced in the pocket. 20.000 switch cycles guaranteed.
- ∴ **Replaceable batteries.** Size AA, or optional rechargeable battery with mini NT table charger\*.



[01]

|   |                     |
|---|---------------------|
| <b>mini3000 Combi Lamp</b>  | <b>2.5V</b>         |
| on mini3000 Battery Handle with 5 disposable blades and batteries | <b>D-001.76.120</b> |
| <b>without handle</b>   | <b>D-001.76.101</b> |
| XHL Xenon Halogen spare bulb                                      | <b>X-001.88.107</b> |
| Box of 100 disposable tongue-blades [01]                          | <b>B-000.12.304</b> |



\*To upgrade to a rechargeable instrument with mini NT table charger, please see page 146.



HEINE mini3000 pocket instruments – a complete range of quality diagnostic instruments for your pocket. Available in black or blue. Please specify the colour you prefer when ordering. If not otherwise specified, we will supply black instruments.