

# Examination Lights

## EXCEPTIONALLY BRIGHT LIGHT FOR EVERYDAY USE

Whether it is the innovative LED Examination light or the tried and true with Halogen, HEINE Examination Lights incorporate the latest and most reliable in illumination technology. They have in common their flexibility and their exceptional brightness.

The three different mounting options by wheeled stand, wall- or clamp-mounting provide ideal solutions for the daily practice environment.



LED NOW IN HEINE QUALITY.





## HEINE EL 10 LED® Examination Light

### Examination light with LED in HQ



The EL10 LED examination light distinguishes itself through an absolutely homogeneous field of light. The light intensity can be regulated continuously between 6,500 and 45,000 Lux. The spot size can be adjusted in a diameter between 8.5 and 14 cm to suit any situation. Both at a working distance of 30 cm.

Controls with ergonomic grip allow for intuitive and fast adjustment of brightness and spot size. The compact illumination head (Ø approx. 60 mm) allows an almost coaxial illumination, particularly in difficult application situations.

- ∴ **Bright and homogeneous.** Perfect illumination in every examination situation.
- ∴ **LED in HQ.** 10 Watt high performance LED with authentic colour rendering (Colour temperature typ. 3,000 Kelvin, CRI > 85).
- ∴ **3 different mounting options.** Wall mounting, clamp for table-top mounting or on a wheeled stand.
- ∴ **Solid construction.** Reliable operation safety and illumination power for many years.
- ∴ **Hygiene.** Easy and efficient cleaning and disinfection.
- ∴ **Coiled goose neck.** Easy and intuitive adjustment.
- ∴ **An electrical control system** monitors all functions and guarantees safety.

- ∴ Continuous adjustment
- ∴ One-hand operation
- ∴ Compact design



**LED HQ**  
LED NOW IN HEINE QUALITY.

EL 10 LED Examination Light	LED
with wall mount [01]	J-008.27.001
with clamp for table-top mounting [02]	J-008.27.002
with wheeled stand, metal base [03]	J-008.27.003
Clamp for table-top mounting for EL10 LED	J-000.27.013
Wheeled stand, metal base with mount for EL10 LED	Y-096.50.010
Mount for wheeled stand	J-000.27.014



[01]



[02]



[03]

## HEINE® HL 5000 Examination Light

### Powerful Universal Halogen Light



With an exceptionally-high light output of up to 120,000 lux at a distance of 30 cm the HL 5000 Examination Light is ideal for examining and recognizing small details. The flexible goose neck construction provides a working radius of 120 cm. Three versions are available: wheeled stand, wall mounting, multi-purpose clamp.

- ∴ **Long-lasting 12V/50W halogen bulb gives up to 120,000 Lux light output at 30 cm.** Very bright, white light.
- ∴ **Extremely compact illumination head.** No obstruction of view, ease of handling.
- ∴ **3 different mounting options: wall, table or wheeled stand.** Versatile.
- ∴ **Halogen bulb lasts at least 4000 hours.** Economical.
- ∴ **Additional UV filter.** Increased safety during examination.
- ∴ **Solid construction with high-specification 120 cm radius goose neck design.** Super-wide radius of operation, easy to position, light stays put in any selected position.
- ∴ **Additional IR-Filter reduces heat.** More comfort for the patient and examiner.

- ∴ Two flexible sections
- ∴ Exceptionally-high light output
- ∴ Up to 120,000 Lux at 30 cm



HL 5000 Examination Light	12V/50W
with wall mount [01]	J-005.27.102
with clamp for table-top mounting [02]	J-005.27.104
with wheeled stand, metal base [03]	J-005.27.106
with lightweight wheeled stand, polycarbonat base [04]	J-005.27.108
Halogen spare bulb 12V/50W [05]	J-005.27.111



[01]



[02]



[03]



[04]



[05]



## HEINE® HL 1200 Examination Light



- ⚡ Two flexible sections
- ⚡ Ultra-slim head
- ⚡ Up to 90,000 Lux at 30 cm



- ⚡ **Ultra-slim head.** More freedom of movement and an unobstructed view for the examiner.
- ⚡ **Built-in temperature filter.** Avoids patient discomfort.
- ⚡ **3 different mounting options: wall, table or wheeled stand.** Versatile.
- ⚡ **Exceeds all relevant quality standards.** Guaranteed protection for the patient even if e.g. a bulb fails.
- ⚡ **“Stay-put” construction of flexible section.** Adjust the light to virtually any position and it “stays-put”, will not sag or droop.
- ⚡ **Long-life 12V/20W lamp.** 4000 hours bulb life.
- ⚡ **Up to 90,000 Lux light output at 30 cm.** Bright, white light.

HL 1200 Examination Light, 12V/20W	black	white
with wall-mounting clamp, can also be fitted to wheeled stand [01]	J-005.27.052	J-005.27.054
with clamp for table-top mounting [02]	J-005.27.056	J-005.27.058
with wheeled stand, metal base [03]	J-005.27.068	J-005.27.067
with lightweight wheeled stand, polycarbonate base [04]	J-005.27.078	J-005.27.077

Halogen spare bulb 12V/20W [05]	J-005.27.075
Wheeled stand (metal base) with mount for HL 1200 Light	Y-096.50.004
Wheeled stand (polycarbonate base) with mount for HL 1200 Light	Y-096.50.009
Mount for HL 1200 Light for wheeled stand	J-000.27.069
Black clamp, only	J-000.27.010
White clamp, only	J-000.27.012



[01]



[02]



[03]



[04]



[05]

230V is the standard mains voltage. If a different voltage is required, e. g. 120V, please specify when ordering.