



Non-Conductive Electrosurgical Instruments & Specula

- 1.** Graves Speculum with smoke tube
Item n°130-100 (Lg.)
Item n°130-101 (Med.)
Item n°130-102 (Sml.)
- 2.** Pederson Speculum with smoke tube
Item n°130-103 (Lg.)
Item n°130-104 (Med.)
Item n°130-105 (sml.)
- 3.** See-More Speculum with smoke tube Item n°130-110 (Lg.)
Item n°130-111 (Med.)
Item n°130-109 (Sml.)
- 4.** Mini See-More Speculum with smoke tube
Item n°130-113
Blade Size: L 70 x W 30mm
- 5.** Kogan Endocervical Speculum
Item n°130-130 (5mm, tips)
Item n°130-131 (4mm, tips)
Item n°133-132 (6mm, tips)
- 6.** Emmett Hook Item n°130-151
- 7.** Iris Hook Item n°130-150
- 8.** Two-Prong Hook Item n°130-154
- 9.** Three-Prong Hook Item n°130-155
- 10.** See-More Speculum with smoke tube
Item n°130-112 (XL)
- 11.** Lateral Vaginal Retractor
Item n°130-141
- 12.** Schroeder Tenaculum
Item N°130-160 (Exposed tips)
Item n°130-161 (coated tips)
- 13.** Collin Speculum with smoke tube
Item n°130-120 (Sml.)
Item n°130-121 (Med.)
Item n°130-122 (Lg.)

- All specula incorporate a smoke evacuation tube.
- All instruments are protected by the highest quality, fully autoclavable blue coating.

All these coated instruments are also available in black (anti-reflexion) for use with lasers.