

PretzelFlex™

part of the EndoFlex™ range


Surgical Retraction System

PretzelFlex™ is the world's first reusable pretzel shaped organ and tissue retractor for use in minimally invasive surgery. It provides improved access and visualisation of surgical sites.

The device has a slim line shaft that can be inserted through a 5mm port. Then, once inside the abdomen, the PretzelFlex™ forms into the innovative pretzel shape.

Key Strengths of PretzelFlex™

- + Impressive strength and stability to support large organs and tissue, ideal for obesity surgery
- + Easily inserted into YelloPort+plus™ 5mm laparoscopic port access system
- + Designed and precision-machined 'in house' to maintain rigidity and shape
- + Smooth slim-line shaft designed to ensure tissue is not penetrated or damaged
- + Innovative segment technology ensures minimal space is required inside the abdomen to form into pretzel shape
- + Ideally suited to be used with FastClamp™ endoscopic clamping system



For further detail about PretzelFlex™ please contact:

sales@surginno.co.uk

+44(0)113 230 7597